CLINICAL IMMUNODIAGNOSTIC AND RESEARCH LAB MACC Fund Research Center, Rm 5072 Phone: 414-456-4165

Fax: 414-456-6487

A Jeffrey Modell in Diagnostic Center for Primary Immunodeficiency

SPECIMEN TRANSPORT GUIDE **COLLECTION GUIDELINES** TEST COLLECTION VOLUME (mL) **TEST SCHEDULE** TURNAROUND **TEST NAME** CODE TUBE REQUIRED MINIMUM DAY APPROX. (days) ABSOLUTE T4 AT4 NaHeparin (green top) 4 Mon - Fri 9AM & 2PM 2 to 3 AUTOIMMUNE LYMPHOPROLIFERATIVE 10 2 to 3 AILYMP NaHeparin (green top) 4 Mon.- Fri. 9AM & 2PM NaHeparin (green top) 10 BRUTON'S TYROSINE KINASE BTK 4 Mon.- Fri. 9AM & 2PM 2 to 3 NaHeparin (green top) COMMON VARIABLE IMMUNODEFICIENCY CVID 10 4 9AM & 2PM 2 to 3 Mon.- Fri CYTOTOXICITY/APOPTOSIS CYTAPO NaHeparin (green top) 15 8 9AM & 2PM 2 to 3 Mon.- Fri 2 to 3 10 HIGM NaHeparin (green top) 4 HYPER IGM Mon.- Fri. 9AM only 8 MENDELIAN SUSCEPTIBILITY TO MYCOBACTE MSMD NaHeparin (green top) 15 Mon., Thurs., Fri. 2PM only 5 to 6 NEUTROPHIL OXIDATIVE BURST **NEUOXB** NaHeparin (green top) 4 1 9AM & 2PM 2 to 3 Mon.- Fri. NaHeparin (green top) NEUTROPHIL PHAGOCYTOSIS 4 PHAGO 1 Mon.- Fri. 9AM & 2PM 2 to 3 PERFORIN GRANZYME PERGRA NaHeparin (green top) 10 4 Mon.- Fri. 9AM & 2PM 2 to 3 10 4 2 to 3 PRIMARY IMMUNODEFICIENCY 1 PID1 NaHeparin (green top) Mon.- Fri. 9AM & 2PM 10 4 9AM & 2PM 2 to 3 PRIMARY IMMUNODEFICIENCY 2 PID2 NaHeparin (green top) Mon.- Fri SEVERE COMBINED IMMUNODEFICIENCY SCID NaHeparin (green top) 2 1 Mon.- Fri. 9AM & 2PM 2 to 3 8 T CELL ACTIVATION TCACT NaHeparin (green top) 15 Mon., Thurs., Fri. 9AM only 5 to 6 T CELL INTERLEUKIN PROLIFERATION TINTL NaHeparin (green top) 15 8 Mon., Thurs., Fri. 2PM only 5 to 6 T CELL MITOGEN PROLIFERATION TMITO NaHeparin (green top) 15 8 Mon., Thurs., Fri. 5 to 6 2PM only 10 4 5 to 6 T HELPER IL17 THIL17 NaHeparin (green top) Mon.- Thurs. 2PM only TOLL - LIKE RECEPTOR TLREC NaHeparin (green top) 15 8 Mon.- Thurs. 2PM only 5 to 6 Mon.- Thurs. 2PM only T REGULATORY - FOXP3 **TREG** NaHeparin (green top) 10 4 5 to 6 NaHeparin (green top) WISKOTT-ALDRICH SYNDROME PROTEIN WASP 10 2 to 3 2 to 3 4 Mon.- Fri. 9AM & 2PM NaHeparin (green top)

SPECIMEN COLLECTION/SUBMISSION/DOCUMENTATION/REQUISITION FORM

10

Mon.- Fri.

9AM & 2PM

- 1 Specimens must be received by the Clinical Immunodiagnostic and Research Lab within 24 hours from the time of blood draw. Recommend timely blood collections and review of the above collection guidelines(see table above).
- 2 Specimens are to be peripheral blood.

X-LINKED LYMPHOPROLIFERATIVE SYNDROMIXLP

- 3 Proper labeling of specimen tube:
 - ▶ Patient's First and Last Name (correctly spelled)
 - ► Medical Record # (Unique Identifier)
 - ▶ Date and Time Drawn
 - Specimen Type or Source
- 4 Identical information must be provided between the specimen collection tube(s) and the submitted requisition form.
- 5 The lab does not bill third parties; therefore, please be sure to provide complete institutional billing information.
- 6 Each specimen must be provided with a COMPLETE Clinical Immunodiagnostic and Research Laboratory requisition form.

IMPORTANT! Requistion is complete when criteria below are me SPECIMEN RELATED INSTITUTION CONTACT ► Name of Ordering Physician Patient Name (First and Last) Medical Record # (Unique Identifier) Direct Contact Name and Date of Birth Phone # for Financial Inquiry Date and Time Drawn • Billing Address Fax # to Send Final Reports Specimen Type or Source Select Test(s) Requested

SHIPPING AND HANDLING

1 Specimens are to be shipped Sunday through Thursday for first priority overnight earliest delivery by 8AM.

Highly Recommend Federal Express: FIRST PRIORITY OVERNIGHT.

Highly Recommend UPS: NEXT DAY AIR EARLY A.M.

- 2 Call the Clinical Immunodiagnostic and Research Lab at 414-456-4165 with tracking # and name of delivery service.
- 3 DOT and IATA requirements for shipping diagnostic specimens:
 - ▶ Place specimen in a leakproof primary receptacle.
 - ▶ Place primary receptacle in leakproof secondary vessel.
 - ▶ Place absorbent material between the primary and secondary packaging.
 - ▶ Primary receptacle or the secondary packaging must meet the required pressure differential of 95 kPa without leaking.
 - ▶ An itemized list of contents must be enclosed between the secondary packaging and the outer packaging.
 - ▶ Package and Air Waybill must show the text: "Diagnostic Specimen" and labeled with the UN3373 designation.
- 4 Specimens must maintain room temperature during transport.

Recommend the use of insulated shipping containers.

5 All delivery carriers (local and non local), deliver to the Medical College of Wisconsin LOADING DOCK.

LOADING DOCK entrance is located off of Watertown Plank Road between 87th & 92nd streets.

Medical College of Wisconsin LOADING DOCK hours are 8AM to 3PM (Monday thru Friday); closed on weekends.

SPECIMEN DELIVERY ADDRESS

Medical College of Wisconsin

Clinical Immunodiagnostic and Research Laboratory

MACC Fund Research Center, Room 5072

8701 Watertown Plank Road

Milwaukee, WI 53226

Effective: 11/2010